TricoCheck® CTRL CONTROLS

DRUGS OF ABUSE, MEDICAMENTS AND EtG IN HAIR

Samples of human hair incorporated and assayed per analyte

CTRL NEG multiCONTROL **CTRL H 20 multiCONTROL** **CTRL EtG CONTROL**
REF SSCTR001270 REF SSCTR002270 REF SSCTR003270

TricoCheck® Screening CTRL-T CONTROLS

DRUGS OF ABUSE AND MEDICAMENTS IN HAIR

Samples of human hair incorporated and assayed per substance class

CTRL-T NEG CONTROL **CTRL-T POS CONTROL** **CTRL-T AMF CONTROL** **CTRL-T THC CONTROL** **CTRL-T BUP CONTROL**
REF SSCTR001270 REF SSCTR002270 REF SSCTR003270 REF SSCTR004270 REF SSCTR005270

CTRL-T BZD CONTROL **CTRL-T KTZ CONTROL** **CTRL-T NIC CONTROL**
REF SSCTR006270 REF SSCTR007270 REF SSCTR008270

TricoCheck® CTRL CONTROLS

DRUGS OF ABUSE, MEDICAMENTS AND ETHYLGLUCURONIDE IN HAIR

EXPECTED USE

The controls TricoCheck® DRUGS OF ABUSE, MEDICAMENTS AND EtG IN HAIR are carried out by samples of human hair, washed and chopped, incorporated with the substances listed inside the sheet SPECIFIC DATA AND VALUES, with the indicated concentrations, to be used in the quality control of dosing procedures of abuse drugs, medicaments and alcohol abuse marker (EtG) in the hair with analytical LC/MS/MS and/or GC/MS and/or screening techniques.

The control TricoCheck® multiCONTROL CTRL NEG is also a negative control for the substances indicated in the sheet SPECIFIC DATA AND VALUES. The TricoCheck® DRUGS OF ABUSE, MEDICAMENTS AND EtG IN HAIR controls can be washed by using normal decontamination, without altering its characteristics relating to the extracting procedure.

TricoCheck® Screening CTRL-T CONTROLS

DRUGS OF ABUSE AND MEDICAMENTS IN HAIR

EXPECTED USE

The controls TricoCheck® Screening CTRL-T DRUGS OF ABUSE AND MEDICAMENTS IN HAIR are carried out by samples of human hair, washed and chopped, incorporated with substances belonging to the classes listed inside the sheet SPECIFIC DATA AND VALUES, with the indicated concentrations, to be used for the check of the whole dosing method of abuse drugs and medicaments in the hair with screening techniques, mainly immunometric type. **The control TricoCheck® Screening CTRL-T NEG is also a negative control for the substance classes listed in the sheet SPECIFIC DATA AND VALUES.**

The controls TricoCheck® CTRL-T DRUGS OF ABUSE AND MEDICAMENTS IN HAIR can be washed by using normal decontamination, without altering its characteristics relating to the extracting procedure.

STORAGE AND STABILITY

The controls TricoCheck® DRUGS OF ABUSE, MEDICAMENTS AND EtG IN HAIR and TricoCheck® Screening DRUGS OF ABUSE AND MEDICAMENTS IN HAIR are stored at room temperature, away from the light and are stable up to the expiration date. They must not be frozen or heated at temperatures above 40°C. Bad storage can affect the proper functioning of the control.

PACKAGE CONTENTS

1 x 270 mg TricoCheck® DRUGS OF ABUSE, MEDICAMENTS AND EtG IN HAIR / TricoCheck® Screening DRUGS OF ABUSE AND MEDICAMENTS IN HAIR. Authentic hair sample, incorporated with the indicated substances, already washed and chopped, ready to be weighed (or pulverized), extracted and analysed with appropriate techniques and extraction and dosing methods. The controls TricoCheck® multiCONTROL CTRL NEG and TricoCheck® Screening Control CTRL-T NEG are negative controls for the substances or the substance classes indicated in the relative sheet SPECIFIC DATA AND VALUES.

CHARACTERISTICS AND PERFORMANCE

By virtue of the particular incorporation technique, the controls TricoCheck® DRUGS OF ABUSE, MEDICAMENTS AD EtG IN HAIR / TricoCheck® Screening DRUGS OF ABUSE AND MEDICAMENTS IN HAIR are particularly homogeneous and stables, they can be further washed – by using normal decontamination techniques – without altering the chemical and physical characteristics and the expected results of the extracting procedure and the following dosing of the incorporated substances. The substances incorporated in the controls TricoCheck® DRUGS OF ABUSE, MEDICAMENTS AND EtG IN HAIR / TricoCheck® Screening DRUGS OF ABUSE AND MEDICAMENTS IN HAIR, can be extracted through techniques described in literature.

INCORPORATED ANALYTES, CONCENTRATIONS, ASSIGNED VALUES, ACCEPTABILITY RANGE

This information is contained in the sheet "SPECIFIC DATA AND VALUES" to be found in the control packing, or it is available on the website: www.comedical.biz/doi/tricocheck - select **REF** and **LOT**. Make sure that the identification information indicated in the sheet "SPECIFIC DATA AND VALUES", matches the corresponding values indicated on the bottle of the control. Any inconsistency can be checked on the website mentioned above.

PRECAUTIONS AND WARNINGS

No study indicating transmission of infections through keratin matrix was found in literature. Nevertheless, since it is not possible to totally exclude the danger of infection, it is necessary to handle the TricoCheck® controls by taking the same precautions used for normal biological samples (use suitable individual protection devices, work with disposable gloves, under fume hood, etc.). Do not use the controls if they are expired. The concentration values of the different analytes or of the substances incorporated in the controls TricoCheck® DRUGS OF ABUSE, MEDICAMENTS AND EtG IN HAIR or TricoCheck® Screening DRUGS OF ABUSE AND MEDICAMENTS IN HAIR may vary between the different lots: always read, before use, the documentation attached to the product.

PROCEDURE

The controls TricoCheck® DRUGS OF ABUSE, MEDICAMENTS AND EtG IN HAIR and TricoCheck® Screening DRUGS OF ABUSE AND MEDICAMENTS IN HAIR can be treated in the same way as a hair sample taken from a human subject (routine) as such they can be washed, even if it is not necessary, to be externally decontaminated.

Weigh the amount (mg) of the control TricoCheck® DRUGS OF ABUSE, MEDICAMENTS AND EtG IN HAIR or TricoCheck® Screening DRUGS OF ABUSE, AND MEDICAMENTS IN HAIR provided by the analytical procedures to be used and apply the appropriate extracting and analytical methods. After analyzing the control, compare the obtained results with the values indicated in the table "SPECIFIC DATA AND VALUES" to be found in the control packing, or on the website: www.comedical.biz/doi/tricocheck - select **REF** and **LOT**.

WARNING

The adulteration of the control, the lack of an adequate preparation or any other non-compliance with label instructions can alter / compromise the performance of the product and invalidate the contents of the label itself explicitly or implicitly. Before opening the bottle, wait for it and its contents to reach room temperature.

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"Buprenorphine detection in hair samples by immunometric screening test: Preliminary experience" - Forensic Science International - F. Svaizer, A. Lotti, M. Gottardi, M.P. Miozzo (Trento, Italy).

"Assessment of Unsuspected Exposure to Drugs of Abuse in Children from a Mediterranean City by Hair Testing" - S. Pichini, O.G. Algar, F. Svaizer, M. Pellegrini, A. Alvarez, M.C. Rotolo, M. Gottardi, R. Pacifici, E. Marchei (Roma, Italy).

"Identification and quantification of 11-nor-D9 tetrahydrocannabinol-9-carboxylic acid glucuronide (THC-COOH-glu) in hair by ultra-performance liquid chromatography tandem mass spectrometry as a potential hair biomarker of cannabis use" Forensic Science International - S. Pichini, M. Pellegrini, R. Pacifici, E. Marchei, F. Svaizer, S. Martello, A. Lotti, M. Gottardi, M. Chiarotti (Roma, Italy).

EXPLANATION OF SYMBOLS

 REF	Catalogue number		Consult instructions for use
 LOT	Batch code		Pay attention
	Use by		Control
	Storage temperature limit/s		CE Mark
	Keep away from light		In Vitro Diagnostic Medical Device
	Dispose of containers and residues suitably		Manufacturer

DISTRIBUTION AND TECHNICAL SUPPORT

Comedical®

Comedical s.a.s. via della Cooperazione, 29 - 38123 Trento (TN) - ITALY
Ph. +39 0461 945876 - Fax+39 0461 944570 - info@comedical.biz - www.comedical.biz

TricoCheck® CTRL CONTROLS

DRUGS OF ABUSE, MEDICAMENTS AND EtG IN HAIR

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“Buprenorphine detection in hair samples by immunometric screening test: Preliminary experience” - Forensic Science International - F. Svaizer, A. Lotti, M. Gottardi, M.P. Miozzo (Trento, Italy).

“Assessment of Unsuspected Exposure to Drugs of Abuse in Children from a Mediterranean City by Hair Testing” - S. Pichini, O.G. Algar, F. Svaizer, M. Pellegrini, A. Alvarez, M.C. Rotolo, M. Gottardi, R. Pacifici, E. Marchei (Roma, Italy).

“Identification and quantification of 11-nor-D9 tetrahydrocannabinol-9-carboxylic acid glucuronide (THC-COOH-glu) in hair by ultra-performance liquid chromatography tandem mass spectrometry as a potential hair biomarker of cannabis use” Forensic Science International - S. Pichini, M. Pellegrini, R. Pacifici, E. Marchei, F. Svaizer, S. Martello, A. Lotti, M. Gottardi, M. Chiarotti (Roma, Italy).

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